## **CSF β2 Transferrin and β Trace Protein**

Accreditation Status:	UKAS Schedule of Accreditation
Date Scheme started:	2011
Clinical Applicability:	The diagnosis of cerebrospinal fluid (CSF) rhinorrhea or otorrhea (leakage of CSF into the nose or ear canal, usually as a result of head trauma, tumour, congenital malformation, or surgery) is often difficult to confirm. CSF B2 Transferrin testing is used to determine the presence or absence of CSF (in serum) in such cases. Beta 2 Transferrin is only found in CSF, ocular fluids and perilymph, therefore it can be used as a marker to determine the presence of CSF in various secretions (typically from the nose and ear)
Analytes:	CSF β2 Transferrin, β Trace Protein <b>(SER/046)</b>
Units for Reporting:	Qualitative: Positive /Negative Quantitative: mg/L
Samples Distributed:	Normal and pathological human serum. Serum based or CSF samples
Number of Distributions per year:	6
Number of Samples per Distribution:	2
Frequency of Distributions:	Every two months as outlined in the Distribution Schedule
Schedule of Analysis:	<b>Data entry</b> is via the web for the submission of results. Data analysis is commenced 28 days after sample dispatch. Late returns are accepted and will contribute to the laboratory's cumulative performance statistics
Schedule of Analysis: Data Analysis:	is commenced 28 days after sample dispatch. Late returns are accepted and will
	is commenced 28 days after sample dispatch. Late returns are accepted and will contribute to the laboratory's cumulative performance statistics Qualitative responses are assessed by MI scoring in relation to the
Data Analysis:	is commenced 28 days after sample dispatch. Late returns are accepted and will contribute to the laboratory's cumulative performance statistics Qualitative responses are assessed by MI scoring in relation to the designated response
Data Analysis: Performance Scoring:	<ul> <li>is commenced 28 days after sample dispatch. Late returns are accepted and will contribute to the laboratory's cumulative performance statistics</li> <li>Qualitative responses are assessed by MI scoring in relation to the designated response</li> <li>MI scoring</li> <li>Laboratory performance is assessed over a running analytical window of 6</li> </ul>